
Conditions for Hybridization of Oligonucleotide Probes

When using oligonucleotides as probes, the aim is to find conditions that are stringent enough to guarantee specificity and sufficiently flexible to allow formation of stable hybrids at an acceptable rate. For DNA molecules more than 200 nucleotides in length, hybridization is usually carried out at 15–25°C below the calculated melting temperature (T_m) of a perfect hybrid. However, as the length of the probe is decreased, the T_m is lowered to the point where it is often impractical to carry out hybridization at $T_m - 25^\circ\text{C}$. Typically, therefore, hybridization with synthetic oligonucleotides is carried out under conditions that are only 5–10°C below the T_m . Although such stringent conditions reduce the number of mismatched clones that are detected with short oligonucleotide probes, they have the less desirable consequence of reducing the rate at which perfect hybrids form.

Hybrids formed between DNA molecules more than 200 nucleotides in length are completely stable for all practical purposes. The chances that such a long stretch of double helix will unwind at temperatures 15–25°C below the T_m are extremely small. However, hybrids (even perfect hybrids) formed between short oligonucleotides and their target sequences at 5–10°C below the T_m are far easier to unwind, and hybridization reactions of this type can be regarded as reversible. This has important practical consequences. Whereas hybrids formed between longer DNA molecules are essentially stable under the conditions used for posthybridization washing, hybrids (even perfect hybrids) involving short oligonucleotides are not. Posthybridization washing of such hybrids must therefore be carried out rapidly so that the probe does not dissociate from its target sequence. For this reason, hybridizations with short oligonucleotides should be carried out under stringent conditions (5–10°C below the T_m) using high concentrations (0.1–1.0 pmole/ml) of probe. When only one or a small number of oligonucleotides (<8) are used as probes, the annealing reaction rapidly reaches equilibrium, and hybridization should therefore be terminated after 3 or 4 hours. More complex mixtures, in which the concentration of each oligonucleotide is comparatively low, require hybridization to be carried out for proportionately longer periods. For example, mixtures of 32 or more oligonucleotides should be hybridized for 1–2 days. Posthybridization washing should be carried out for brief periods of time, initially under conditions of low stringency and then under conditions of stringency equal to those used for hybridization.

CALCULATING MELTING TEMPERATURES FOR PERFECTLY MATCHED HYBRIDS BETWEEN OLIGONUCLEOTIDES AND THEIR TARGET SEQUENCES

When using single oligonucleotides that match the target sequence perfectly, hybridization conditions can easily be derived from the calculated T_m of the hybrid. For oligonucleotides shorter than 18 nucleotides, the T_m of the hybrid can be estimated by multiplying the number of A + T residues in the hybrid by 2°C and the number of G + C residues by 4°C and adding the two numbers (Itakura et al. 1984). However, this method overestimates the T_m of hybrids involving longer oligonucleotides.

A different approach has been taken by E. Fritsch (unpubl.), who found that the equation originally used to calculate the relationship between G + C content, ionic strength of the hybridization solution, and the T_m of long DNA molecules (Bolton and McCarthy 1962):

$$T_m = 81.5 + 16.6(\log_{10}[\text{Na}^+]) + 0.41(\text{fraction G} + \text{C}) - (600/N),$$

where N = chain length, predicts reasonably well the T_m for oligonucleotides as long as 60–70 nucleotides and as short as 14 nucleotides.

This formula only works for Na^+ concentrations of 1 M or less.

ESTIMATING THE EFFECTS OF MISMATCHES

Perhaps surprisingly, the classic formula (Bonner et al. 1973) to calculate the effect of mismatches on the stability of long DNA hybrids holds reasonably well for hybrids involving short oligonucleotides: For every 1% of mismatching of bases in a double-stranded DNA, there is a reduction of T_m by 1–1.5°C. However, the precise effect of mismatches depends on the G + C content of the oligonucleotide and, even more critically, on the distribution of mismatched bases in the double-stranded DNA. Mismatches in the middle of the oligonucleotide are far more deleterious than mismatches at the ends. Therefore, the method of estimation given above should only be used as a rough guide until a systematic study of all types of mismatches in a variety of contexts leads to more precise methods of estimation. If appropriate target DNA has been cloned, the effect of mismatches on T_m can be determined empirically (see pages 11.55–11.57).

HYBRIDIZATION OF POOLS OF OLIGONUCLEOTIDES

It is easy to calculate accurately the T_m of a perfectly matched hybrid formed between a single oligonucleotide and its target sequence. However, when using pools of oligonucleotides whose members have greatly different contents of G + C, it is impossible to estimate a consensus T_m . Because it is not possible to know which member of the pool will match the target sequence perfectly, conditions must be used that allow the oligonucleotide with the lowest content of G + C to hybridize efficiently. Usually, conditions are chosen to be 2°C below the calculated T_m of the most A/T-rich member of the pool (Suggs et al. 1981b). However, the use of such "lowest common denominator" conditions can lead to a number of false positives, because mismatched hybrids formed by oligonucleotides of higher G + C content may be more stable than a perfectly matched hybrid formed by the correct oligonucleotide. In most cases, this problem is not serious, since the number of positive clones obtained by screening cDNA libraries with pools of oligonucleotides is usually quite manageable. It is therefore possible to easily distinguish false positives from true positives by another test (e.g., DNA sequencing or hybridization with a second pool of oligonucleotides corresponding to another segment of amino acid sequence).

In those cases when the number of positives is unacceptably high, it may be worthwhile to consider using hybridization solvents that contain the quaternary alkylammonium salts tetraethylammonium chloride (TEACl) or tetramethylammonium chloride (TMACl) instead of sodium chloride (Melchior and von Hippel 1973; Jacobs et al. 1985, 1988; Wood et al. 1985; Gitschier et al. 1986; Wozney 1989). In these solvents, the T_m of a hybrid is independent of its base composition and dependent primarily on its length. Thus, by choosing a temperature for hybridization appropriate for the lengths of the oligonucleotides in a pool, the effects of potential mismatches can be minimized.

It is important to obtain an accurate estimate of the T_m in TMACl or TEACl before using pools of oligonucleotides to screen cDNA or genomic DNA libraries. Jacobs et al. (1988) measured the T_i (the irreversible melting temperature of the hybrid formed between the probe and its target sequence) as a function of chain length for a number of oligonucleotides of differing G + C content in solvents containing either sodium or tetramethylammonium ions. Hybrids involving oligonucleotides 16 and 19 nucleotides in length melt over a smaller range of temperature in solvents containing TMACl than in solvents containing sodium salts (3°C for TMACl vs. 17°C for SSC when hybridizing 16-mers; 5°C for TMACl vs. 20°C for SSC when hybridizing 19-mers). For 14-mers, the effect is much less dramatic (7°C for TMACl vs. 9°C for SSC). Similar, but less extensive, data are available for solvents containing TEACl (Jacobs et al. 1988).

The optimal temperature for hybridization is usually chosen to be 5°C below the T_i for the given chain length. The recommended hybridization temperature for 17-mers in 3 M TMACl is 48–50°C; for 19-mers, it is 55–57°C; and for 20-mers, it is 58–66°C. Three points are worth emphasizing. First, the T_i s of hybrids are uniformly 15–20°C higher in solvents containing TMACl than in solvents containing TEACl. The higher T_i in solvents containing TMACl allows hybridization to be performed at temperatures that

suppress nonspecific adsorption of the probe to solid supports (such as nylon membranes), resulting in lower nonspecific backgrounds. Second, hybridization solvents containing TMACl do not have significant advantages over those containing sodium ions until the length of the oligonucleotide exceeds 16 nucleotides. Finally, the data have been extensively examined for 16-mers, 19-mers, and, in previous studies, for much longer DNA molecules (Melchior and von Hippel 1973). It is currently an untested assumption that the same beneficial effect will be seen for DNA molecules of all intermediate lengths.